



Manual

Terahertz generator

Terastim Headband ProTM

The Terastim series is a natural way to activate the body's own resources

Människans Resurser AB

Terastim Headband PRO with magnets

Brain activator

This headband works in the same way as the regular headband with a built-in chip that activates the production of nitric oxide (NO) in the blood vessel walls and makes them more elastic. The production of NO leads to an improved blood circulation.

In this headband, in addition to the terahertz-producing chip, there are also weak magnets that enhance the effect of the chip, which speeds up the improvement of the blood circulation. This headband is particularly effective for acute situations such as headaches, fatigue and even body aches thanks to improved control abilities of the brain. In addition, this headband with magnets is the best health tool during weather changes, something that usually provokes spasms in the blood vessels, which can cause various painful effects.

When the blood circulation in the brain is improved, many other positive effects may occur - such as reduced headaches, normalization of blood pressure, improved memory etc.

HOW TO USE

You use the headband by placing the logo „Terastim” against the forehead - the chip is on the opposite side of the logo, on the back of the head. This improves blood supply to the pituitary and the hypothalamus, which strengthens the body's regulatory mechanisms.



Use the headband with magnets for 15 min. per session, 3-4 times per day. In acute situations, you can use it up to 30 min. per session.

Contraindications: stroke, cerebral infarction.

You can use hand washing for the headband.

Manufacturer's warranty is 1 years.

Contact

Människans Resurser AB



Byängsgränd 14, tr 3,
120 40 Årsta



mail@mr-ab.se
www.terastim.com



08 334908